

COVID-19 Drive-Through Specimen Collection Station Solutions



Challenge

Optimize COVID-19 drive-through specimen collection with safety, speed and accuracy.


Solution

Leverage mobile barcode technology for data collection, positive patient ID, accurate labeling and printing.

Good

Laptop-Based/Electronic Health Record Configuration


Approx 5-6 Minutes | No Physical Contact | Additional Footsteps | Requires Internet Connection



Better

Turn-Key Mobile Configuration


Under 3 Minutes | No Physical Contact | No Manual Entry | No Internet Connection



Plus

Emergency Wristband Kit

For Patients in Critical Need
Can Be Used With Any Configuration



For more information, please visit zebra.com/covid19