

Mother's Milk

Mobile Care from Mother to Baby

- ✓ Improve Infant **Safety**
- ✓ **Low-Cost** Implementation
- ✓ **Reduce the Risk** of Wrong Feedings
- ✓ **Reporting** and **Compliance** Management
- ✓ **Non-Invasive Scanning** at Patient Bedside
- ✓ Process **Improvement**



“ ...Since the implementation of Mother's Milk we have had no breast milk errors. It has brought a higher level of quality to our patient care and I would strongly recommend this application to other hospitals who offer expressed breast milk to their patients. ”

*Deana Kuhlman, RNC-NICU, BSN
Nursing Supervisor
Pomona Valley Hospital Medical Center*

How It Works



Step 1.

Scan baby's wristband into the mobile solution.



Step 2.

Print matching labels and apply to infant milk bottles.



Step 3.

After the bottles have been filled, scan the labels into the mobile solution.



Step 4.

To complete a feeding, scan the baby wristband to verify a match between bottle and baby.



Step 5.

Access to real-time reporting.

Mother's Milk is a mobile solution developed with caring nurses in mind. Its verification system protects both patients and nursing staff by preventing feeding errors from occurring inside the NICU. The Mother's Milk platform integrates into existing Health Information Systems (EMR) and combines patient identification and breast milk verification using the mother-to-infant match functionality. Through secure login verification, administrators can manage users and generate reports for successful, denied, and overridden scanning verifications.



ANTON ANSALMAR / Rapid Healthcare Inc.

C 714.600.2848 / **T** 949.825.7013 / **E** anton@rapid-healthcare.com

W www.rapid-healthcare.com / **A** 5151 California Ave, Suite 100, Irvine, CA 92617